# **MICROLIT**

**Motor Operated Burett** 3 Calibrated Pre-set Spe

# **Microlit India**

629 Pakramau, Kursi Road, Lucknow 226026, India Phone: +91 9918625629, Email: info@microlit.com www.microlit.com

# **Microlit USA**

33 Wood Avenue South, Suite 600, Iselin, NJ 08830, USA Phone: +1 732 321 0852, Email: info-usa@microlit.com www.microlit.us

# **Microlit LATAM**

Business Place Ibirapuera, Alameda dos Jurupis,1005, Indianópolis, São Paulo, Brasil Phone: +55 11 93269-5932, MICROLITE-BURETTE



# EXPERIENCE PRECISION

# TABLE OF CONTENTS

1.	Intended Use of the Instrument
2.	Product Specifications
3.	Safety Instructions
4.	Functions and Limitations of Use
5.	Operating Exclusions
6.	Storage Conditions
7.	List of Recommended Reagents
8.	First Steps
9.	Initializing
10.	Setting Up
11.	Computer Connectivity
12.	Operating the instrument
13.	Data Transfer Setup
14.	Error Limits
15.	Calibration
16.	Other Settings
17.	Cleaning and Product Maintenance
18.	Troubleshooting
19.	Return for repair
20.	Calibration Guidelines
21.	Warranty
22.	Disposal

1			
1			
2			
3			
4			
5			
6			
13			
16			
18			
20			
23			
29			
33			
34			
36			
38			
39			
39			
40			
40			
40			

#### Intended Use of the Instrument

MICROLIT E-BURETTE is a state-of-the-art Motor Operated Burette with sophisticated features and functionality. It is intended for use in laboratories and the field for easy, safe and accurate titrations.

#### **Product Specifications** 2.

Performance Data	10 ml	25 ml	50 ml		
Unit Type	Microprocessor-controlled motorized piston lift system with recirculation valve				
Volume Ranges	0.005 ml - 99.99 ml (Max piston lift 10 ml, automatic addition)	0.01 ml - 250 ml (Max piston lift 25 ml, automatic addition)	0.01 ml - 500 ml (Max piston lift 50 ml, automatic addition)		
Volume Accuracy	A= 0.10%, CV = 0.050%	A = 0.07%, CV = 0.024%	A= 0.05%, CV = 0.020%		
Operating Temperature	+10 °C to +40 °C (from 50 °F to 104 °F)				
Power Supply	18 V, 3 A				
Battery Operation	One complete battery charge gives 3 hours (approx.) of continuous operation, (Battery capacity indicated on TFT display and charging time 6.5 hrs.)				
Quality Control	In compliance with DIN EN ISO 8655-3				
Touchscreen	4" TFT screen				
PC Interface	Mini USB Cable				
Languages	Available in 4 languages (English, Spanish, German and French)				

#### 3. **Safety Instructions**

This instrument may sometimes be used for operations involving hazardous materials and equipment. It is beyond the scope of this manual to address all of the potential risks associated with its use in such applications. It is the sole responsibility of the user to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations, prior to the usage of this instrument.

- 1 Please read the following instructions carefully.
- Read and understand this Operation Manual thoroughly before ٠ using the instrument.
- Follow general and safety instructions for hazard prevention. E.g., always wear ٠ protective clothing, protective gear for the eyes and hands.
- Carefully observe the specifications provided by reagent manufacturers. •
- When dispensing inflammable fluids, avoid the built up of static charge. • Make sure that you do not dispense into plastic vessels and do not wipe the instrument/equipment with a dry cloth.
- Use the instrument only for dispensing liquids, with strict regard to the defined • limitations of use and operating constraints. (Refer to section 4)
- Observe operating exclusions. When in doubt, contact the • manufacturer or supplier.
- Always use the instrument in such a way that neither the user nor any other person is in danger. While dispensing, the discharge tube must always point away from you or any other person. Only dispense into appropriate vessels and avoid splashes.
- Do not press the piston when the discharge tube closure is attached. •
- Do not remove the discharge tube while the dispensing cylinder is being filled. •
- Reagents can accumulate in the cap of the discharge tube. Thus, clean it regularly. •
- Do not carry the mounted instrument by the cylinder sleeve or the valve block. • Breakage or loosening of the cylinder may lead to personal injury.

1

- Use only original accessories and spare parts.
- Do not attempt to make any technical alterations. Do not dismantle • the instrument any further than is described in the Operation Manual.
- Always check the instrument for visual damage before use. •
- If there is any sign of a potential malfunction (e.g. piston difficult to move, stuck valve or leakage), immediately stop dispensing. Consult the 'Troubleshooting' section of this Operation Manual and contact the manufacturer if needed. (Refer to section 18)

#### **Functions and Limitations of Use**

MICROLIT E-BURETTE is designed for dispensing liquids directly from the reservoir bottle. It is calibrated in accordance with the guidelines of the DIN EN ISO 8655 - 3. When the instrument is correctly used, the dispensed liquid comes into contact with only the following chemically resistant materials: PTFE, FEP and Borosilicate glass.

#### CE CE Marking Α.

This sign certifies that the product meets the requirements of the EC directive and has been tested according to the specified test methods.

#### Β. Physical Limits

MICROLIT E-BURETTE is designed for titrating liquids, observing the following physical limits:

- Maintain instrument and reagent between +10 °C to +40 °C (from 50 °F to 104 °F).
- Maintain reagent vapor pressure up to 600 mbar (maximum). • Aspirate slowly above 300 mbar in order to prevent the liquid from boiling.
- Kinematic viscosity upto 500 mm2/s. • (dynamic viscosity [mPas] = kinematic viscosity [mm2 /s] x density [g/cm<sup>3</sup>])

- Relative air humidity between 20% and 90%. ٠
- Use fluids with density upto 2.2 g/cm<sup>3</sup>.
- **Operating Limitations** C.
- Chlorinated and fluorinated hydrocarbons or chemical combinations which • form deposits may make it difficult for the piston to move or jam it. If the piston becomes too difficult to move, immediately clean the instrument. (Refer to section 17)
- While dispensing inflammable fluids, avoid the built up of static charge. • Make sure that you do not dispense into plastic vessels and do not wipe instruments with a dry cloth.
- MICROLIT E-BURETTE is designed for general laboratory applications and • complies with relevant standards, like DIN EN ISO 8655-3. Please check the compatibility of the instrument for a specific application (e.g. trace material analysis, food sector etc.). Approvals for specific applications, like production and administration of food, pharmaceuticals and cosmetics are not available.

#### **Operating Exclusions** 5.

- Do not use the instrument with: •
- Liquids that attack FEP, PFA and PTFE (e.g. dissolved Sodium Azide\*).
- Liquids that attack borosilicate glass (e.g. Hydrofluoric acid). •
- Hydrochloric acid > 40%.
- Tetrahydrofuran or Oxolane. ٠
- Trifluoroacetic acid.
- Explosive liquids (e.g. Carbon Disulfide).
- Suspensions (e.g. of charcoal) as the solid particles may clog or • damage the instrument.

- Liquids that attack PP (cap).
- An aggressive atmosphere, e.g. HCl fumes.

\* The permissible concentration of dissolved Sodium Azide is not more than 0.1%.

#### 6. Storage Conditions

Store the instrument and accessories in a clean, cool and dry place. Storage temperature from -20 °C to +50 °C (from -4 °F to 122 °F) is advised with relative humidity between 5% and 95%.

## 7. List of Recommended Reagents

#### Chemicals from A to Z

The following tables enlist the most frequently used reagents. They provide useful information for the safe and efficient use of MICROLIT E-BURETTE. All the safety precautions and recommendations within this Operation Manual must be followed carefully.

#### Code explanations


- A = Good resistance
- B = Acceptable with limitations
- C = Not recommended
- 1 = Acid vapours (better resistance with lower concentration). Do not leave the instrument on the bottle.
- 2 = Risk of damage, softening or discoloration of external parts through vapours. Do not leave the instrument on the bottle.
- 3 = Chemical degradation of glass parts (plunger/barrel).

# List of Reagents

#### Chemicals A - Z А Acetaldehyde (Ethanal) А Acetic acid 96% А Acetic acid 100% (glacial) B/2 Acetic anhydride B/2 Acetone (Propanone) B/2 Acetonitrile (MECN) Α Acetophenone B/2 Acetyl Chloride B/2 Acetylacetone А Acrylic acid А Acrylonitrile B/2 Adipic acid А Allyl alcohol А Aluminum chloride А Amino acids А Ammonia 20% B/2 B/2 Ammonia 20-30% Ammonium chloride А Ammonium fluoride Α Ammonium molybdate А Ammonium sulfate А Amyl alcohol (Pentanol) А Amyl chloride (Chloropentane) B/2 Aniline А Ascorbic acid А n-Amyl acetate B/2 В Barium chloride А Benzaldehyde А Benzene B/2 Benzine А Benzoyl chloride B/2 Benzyl alcohol А Benzyl chloride B/2 Bis(2-ethylhexyl) phthalate B/2 Boric acid 10% А Bromine C/2 Bromobenzene B/2 Bromonaphtalene А Butanediol А Butanol А Butyl acetate B/2 Butyl methyl ether B/2 Butylamine B/2 Butyric acid B/2 С Calcium carbonate А Calcium chloride А Calcium hydroxide А Calcium hypochlorite А

# List of Reagents

#### Chemicals A - Z С Carbon disulfide B/2 Carbon tetrachloride B/2 Chlorine dioxide B/2 Chlorine water B/2 Chloro naphthalene B/2 Chloroacetaldehyde 45% А Chloroacetic acid А Chloroacetone B/2 Chlorobenzene B/2 Chlorobutane B/2 B/2 Chloroethanol Chloroform (Trichloromethane) B/2 Nitro-hydrochloric acid (Aqua regia) B/2 Chlorosulfonic acid B/2 Chlorosulfuric acid 100% B/1/2 Chromic acid 100% B/1/2 Chromosulfuric acid 100% C/1/2 Citric acid А Copper fluoride А Copper sulfate А Covi-Ox-T70/ Mixed Tocopherol А Cresol А Cumene (Isopropylbenzene) B/2 Cyanoacrylate А Cyclohexane B/2 Cyclohexanone B/2 Cyclopentane B/2 D 1,2-Diethylbenzene B/2 1,4-Dioxane (Diethylene dioxide) B/2 1-Decanol А Decane А Di-(2-ethylhexyl) peroxydicarbonate B/2 Dibenzyl ether B/2 Dichloroacetic acid А Dichlorobenzene А Dichloroethane А Dichloroethylene B/2 Diesel oil (Heating oil) А Diethanolamine А Diethylamine B/2 Diethylene glycol А Diethylether B/2 Dimethylacetamide А Dimethyl sulfoxide (DMSO) B/2 Dimethylaniline А Dimethylformamide (DMF) B/2



#### List of Reagents Chemicals A - Z E Ethanol Α B/2 Ethanolamine Ether B/2 Ethyl acetate B/2 Ethylbenzene B/2 Ethylene chloride B/2 Ethylene diamine А Ethylene glycol А Fluoroacetic acid B/2 Formaldehyde (Formalin) А Formamide А Formic acid Α Gamma-butyrolactone А Gasoline Glycerin <40% B/2 А Glycolic acid 50% А Heating oil (Diesel oil) А Heptane Hexane А Α Hexanoic acid Hexanol A A Hydriodic acid B/2 А Hydrobromic acid Hydrochloric acid 20% (HCI) Hydrochloric acid 37% (HCI) Hydrofluoric acid (HF) Α B/1 C/3 Hydrogen peroxide А lodine А lodine bromide C/2 C/2 lodine chloride

G

Н

Isoamyl alcohol

Isobutanol

Isooctane

Lactic acid

Kerosene

Methanol

Μ

Isopropanol

Isopropyl ether

Iso-propylamine

2-Methoxyethanol

Methoxybenzene (Anisol) Methyl benzoate

Methyl chloride (Chloromethane)

M	D/0
Methyl ethyl ketone (MEK/Butanone)	B/2
Methyl formate	A
Methyl iodide (lodomethane)	B/2
Methyl methacrylate (MMA)	B/2
Methyl propyl ketone (2-Pentanone)	A
Methyl tert-butyl ether	B/2
Methylene chloride (Dichloromethane) (DCM)	B/2
Methylpentanone	A
Mineral oil (engine oil)	A
Monochloroacetic acid	A
N	
Nitric acid 100%	A
Nitric acid 30-70%	A
Nitric acid dil. <30%	A
Nitrobenzene	B/2
Nitromethane	B/2
N-methyl-2-pyrrolidone (NMP)	A
0	
Octane	A
Octanol	A
Oil (vegetable, animal)	B/2
Oil of turpentine	B/2
Oleic acid	A
Oleum (Fuming Sulfuric acid)	A
Oxalic acid	Α
P	
Pentane	B/2
Peracetic acid	A
Perchloric acid 100%	B/2
Perchloric acid diluted	A
Perchloroethylene	B/2
Petroleum	B/2
Petroleum ether / spirit	B/2
Phenol	Α
Phenylethanol	B/2
Phenylhydrazine	B/2
Phosphoric acid 100%	Α
Phosphoric acid 85%	A
Piperidine	B/2
Potassium chloride	Α
Potassium dichromate	А
Potassium dihydrogen phosphate	Α
Potassium hydroxide	А
Potassium iodide	А
Potassium permanganate (persulfate)	А
Potassium peroxydisulfate	A
Potassium sulfate	А
	A
Propionic acid (Propanoic acid)	A

# List of Reagents

Chemicals A - Z

Α

А

Α

А

B/2

B/2

А

А

А

A B/2


B/2


B/2




Chemicals A - Z	
Propylene oxide	A
Picric acid (Trinitrophenol)	B/2
P	
Pyridine	B/2
Pyruvic acid	А
R	
Resorcin	A
S	
Salicylaldehyde	A
Scintilation fluid	A
Silver acetate	A
Silver nitrate	A
Sodium acetate	A
Sodium chloride (kitchen salt)	A
Sodium dichromate	A
Sodium fluoride	A
Sodium hydroxide 30%	A
Sodium hypochlorite	A
Sodium thiosulfate	A
Sulfonitric acid 100%	B/2
Sulfur dioxide	B/2
Sulfuric acid 100%	B/2
Sulfuric acid <10%	A
Sulfuric acid (10-75%)	B/1
Sulfuric acid (Cold conc.)	A
Sulfuric acid (Hot conc.)	B/2
T	
1,1,2-Trichlortrifluoroethane	B/2
Tartaric acid	A
Tetrachlorethylene	B/2
Tetrahydrofuran (THF)	B/2
Tetramethylammonium hydroxide	A
	B/2
Trichlorethylene	B/2
Trichloroacetic acid	B/2
Trichlorobenzene	B/2 B/2
Trichloroethane Triethanolamine	B/2
Triethylamine Triethylama glycal	A
	A
Trifluoroacetic anhydride (TFAA)	B/2
Trifluoromethane (Fluoroform)	B/2
Urea	A
X	~
∧ Xylene	B/2
Z	DIZ
Zinc chloride 10%	A
Zinc sulfate 10%	A










Adjust the length of the Telescopic tube (FEP tube) to fit your particular reservoir. If you require a longer tube, it is provided on request. (Fig. 10.1) Fix the Telescopic tube as shown. The tube's ends have different diameters - fix it from the wider end. (Fig. 10.2) Choose the correct adapter for the bottle. The threaded platform base of MICROLIT E-BURETTE has a 30 mm screw thread. Five adapters are supplied to suit containers with a 28, 32, 38, 40, 45 mm and 30 mm (inbuilt adapter) screw neck. (Fig. 10.3) Fix the adapter on the reservoir bottle by screwing it in clockwise direction. (Fig. 10.4)

18





Connect the cable from the instrument housing to the port at the back of the control panel with a firm push, as shown. Check for the cable position. The flat surface of the connector with an arrow mark will face down. (Fig. 10.5)



If charging is required, insert the charging pin at the back of the control panel and connect the power cord with the AC adaptor firmly and then plug the two pin power cord to the 220 V power supply. (Fig. 10.6)



The battery indicator will show the charging symbol. (Fig. 10.7)



The E-Burette is now ready to be switched on. Press the ON button (of the ON/OFF button) at the back of the control panel to switch on the control panel. (Fig. 10.8)





When mounted to a reagent bottle, always carry the instrument as shown and always store it in the upright position. (Fig. 10.9)

#### 11. Computer Connectivity

Insert the pen drive provided with the package in your PC/Laptop for driver files and software. Please note that there are two set up files in the pen drive.

- CDM21228
- E-Burette Setup Files

Procedure to install the driver & software:


1. Installing the CDM21228 driver files:



Double click the CDM21228 setup icon. (Fig. 11.1)



Press the Extract tab from the FTDI CDM driver window. (Fig. 11.2)





#### 3.




- Turn the knob to re-circulation mode and click the tick mark ( $\checkmark$ ) on the control panel screen when ready. (Fig. 12.5)
- 4. The E-BURETTE will purge automatically. If the device is still not properly purged, go to step 1 and repeat the procedure till no large air bubble is visible below the piston.

A few air bubbles up to 1 mm in size are permissible.

#### Titration

▲ Wear protective clothing, protective gear for the eyes and hands. Liquid may accumulate in the cap. Follow all safety instructions and observe limitations of use and operating constraints. (Refer to section 4 and 5)



Remove cap from the discharge tube. (Fig. 12.6)



Click on the Titration mode on the home screen. (Fig. 12.7) 3. The screen on the control panel will warn the user to turn the knob to titration mode.



Turn the knob clockwise to set at titration mode and click the tick (✓) mark on the screen of the control panel when ready. (Fig. 12.8)

5. Place the discharge tube orifice against the inner wall of a suitable receiving vessel.





Now fill the barrel from the titration screen by clicking the FILL button (Fig. 12.9 - left part). By pressing the STOP button, the filling volume can be stopped as and when you need to. (Fig. 12.9 - right part)

7. Three dispensing speeds are provided, including drop-wise dispensing, which allows the user to achieve the end point very accurately.



Fast mode In this mode, the dispensing speed is 2.2 ml/sec (Fig. 12.10)





#### Sensitive Media


For light-sensitive media (e.g., lodine, Potassium Permanganate and Silver Nitrate solutions), we recommend the use of 'Amber Coloured Windows' or the inspection windows provided in the box.



To replace the inspection window, unclip the default inspection window by pressing it from its top and then removing it. (Fig. 12.14)



Place the Amber Coloured Window in the slot provided at the bottom of the housing and press it down until you hear a click sound. (Fig. 12.15)



# Memory

1. Click ADD from the control panel to save the readings from the control panel to the internal memory.



To recall the previous recorded reading (s), press the RECALL button. (Fig. 13.8 A)



Please note that the recall grid displays the last 20 recorded readings. The LAST saved reading will be displayed on the number 1 position of the table. (Fig. 13.8 B)

3.



For quick recall of the last reading, press the LAST button. (Fig. 13.9)

Please note that the last saved readings will be displayed till the time LAST button is pressed.

## Creating a File Name



1. Select the port from the port list dropdown menu. 2. Click Open port. 3. The Software captures the Model No., Sr. no. and volume range. 4. After transferring the data from control panel to the software window we can save the file by clicking the "Save as Excel File". The data recorded while operating the E-BURETTE can be printed directly from the software by clicking the PRINT tab. (Fig. 13.11) For this, select the printer and print the readings. You can also access the saved .xls file to print. To start afresh for new readings on the software, click on the REFRESH tab. (Fig. 13.12) The readings on the grid shall be deleted and new readings can now be transferred in it. 32

# Printing the data



Refresh





- go to the home position.
- 3. Press the ON/OFF button to switch OFF the control panel.

## 14. Error Limits

The error limits (Accuracy and Coefficient of Variation) mentioned in Product Specifications (Page 1) are in accordance with the nominal capacity (or maximum volume) indicated on the instrument.

These are obtained by using the instrument with distilled water at equilibrium, ambient temperature of 20 °C, while operating the device smoothly and steadily.

The error limits are well within the limits of DIN EN ISO 8655-3.

Capacity (ml)	Increment	Accuracy ±%	y ±µl	CV ±%	±μΙ
10	5 µl	0.10	10	0.050	5
25	10 µl	0.07	18	0.024	6
50	10 µl	0.05	25	0.020	10

# Calibration

#### User Calibration Procedure

The E-BURETTE has been laboratory calibrated at its nominal volume. However, due to changes in environmental conditions and the viscosity of the media which is dispensed, we recommend gravimetric testing every 3-12 months. Gravimetric volume testing according to DIN EN ISO 8655-6 is performed as follows:



right of the Home Screen.



(Fig. 15.2)



Click the settings icon on the bottom (Fig. 15.1) From the settings menu, choose Calibration. Fill the nominal volume of double distilled de-ionized water by clicking the FILL button on the calibration screen of the operating panel. (Fig. 15.3)

# 4.



Dispense the filled liquid by clicking on the DISPENSE button on the calibration screen of the control panel. (Fig. 15.4)

5. Measure the dispensed liquid on a balance and enter the value in mg using the provided electronic keyboard. Click the tick mark ( ) to proceed.

Note: CAL label will start appearing on top of all screens.

6. Repeat this procedure until the nominal volume is achieved on the electronic balance.

#### Factory Reset



The factory reset is for the MICROLIT factory settings which is password protected and intended to be used at the manufacturing facility. (Fig. 15.5)

If the reset button is pressed, the user calibrated settings shall be reset with the default factory settings & the CAL label will disappear from the screens.

# 16. Other Settings

#### Control panel brightness setting

User can change brightness from setting screen on the control panel.



Click on the settings icon from the home screen. (Fig. 16.1)



In the settings windows, click the BRIGHTNESS icon. (Fig. 16.2)



From the BRIGHTNESS window, select the level of brightness and tick mark to save it. (Fig. 16.3)

# Language Settings

User can select the preferred language from the language settings.



Click on the settings icon from the home screen. (Fig. 16.4)



Select the preferred language from the tab by clicking on it and tick mark ( $\checkmark$ ) to save it.

MICROLIT E-BURETTE supports four languages -English, German, Spanish and French.

(Fig. 16.5)

#### 17.

**Cleaning and Product Maintenance** Whenever cleaning is required, run the burette under distilled water. Fill the barrel completely with distilled water and press the fast dispense button to dispense the water completely. Procedure to Disassemble the Delivery Pipe: Unscrew the chuck nut by rotating it in anticlockwise direction and pull out the delivery pipe. (Fig. 17.1) Clean the delivery pipe with de-ionized water. Procedure to Assemble the Delivery Pipe: First push the delivery pipe into the lower housing till it stops going in further. (Fig. 17.2) Screw the chuck nut to complete the assembly by turning it in clockwise direction. (Fig. 17.3) Note: Use rubbing alcohol on a cloth or cotton to clean the external body.







#### 18. Troubleshooting

Problem	Possible cause	Solution
Piston difficult to move	Formation of crystals or deposition of dirt	Perform a cleaning cycle.
Filling not possible	Filling valve stuck	Clean the filling valve. If the valve ball is stuck, use a 200 µl plastic tip to loosen it. (Fig. 18.1)
Air bubbles in the instrument	Purging not complete Filling tube loose or damaged Filling tube does not dip into the liquid	Purge the instrument again. Fasten the telescopic filling tube firmly. If necessary, cut the tube off by approx. 1 cm from the top or replace it. Fill up the bottle, or correctly adjust the length of the telescopic filling tube.
The volume delivered is smaller than that indicated	The instrument has not been completely primed	Prime the instrument again.
The instrument doesn't indicate any function	Internal error	Perform a restart.

#### 19. **Return for repair**

- 1. Clean and decontaminate the instrument carefully.
- 2. Always mention the exact description of the type of malfunction and the media used. Please note that if information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender. 3.

#### 20. **Calibration Guidelines**

volumetric instruments. We recommend checking the volume every 3-12 months. The interval depends on the specific requirements of the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter.

#### 21. Warranty

the instructions of the operating manual.

#### 22.



Disposal

ISO 9001 and GLP-guidelines require regular examinations of your MICROLIT shall not be liable for the consequences of improper handling, use, servicing, operation or unauthorized repairs of the instrument or the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow We are not liable for damage resulting from any actions not described in the operating manual or non-original spare parts or components being used. The adjoining symbol means that storage batteries and electronic levices must be disposed of separately from household trash 🐨 🚺 (mixed municipal waste) at the end of their service life. According to the Directive 2002/96/EC of the European Parliament and of the Council on Waste Electrical and Electronic Equipment (WEEE) published on 27 January 2003, electronic equipment requires disposal according to the Batteries contain substances that can have harmful effects on the environment and human health. Therefore, according to the Directive 2006/66/EC of the European Parliament and the Council on Waste Batteries of 6 September, 2006, batteries require disposal according to the relevant national disposal regulations. Do not short-circuit the battery to discharge it. Subject to technical modifications

- relevant national disposal regulations.
- Dispose the batteries only when they are completely discharged. without notice. Errors accepted.