

MICROLIT LILPET PRO

In-built Tip Ejector

One-step Plunger Movement

Fully Autoclavable

Fixed Volume

 MICROLIT

SPECIFICATIONS AND ORDERING INFORMATION	Model	Volume	Acc. $\pm\%$	CV $\pm\%$
	LILP-S-5	5 μl	5	3
	LILP-S-10	10 μl	2	1.5
	LILP-S-20	20 μl	2	1
	LILP-S-25	25 μl	2	1
	LILP-S-50	50 μl	2	1
	LILP-S-100	100 μl	1.5	0.6
	LILP-S-200	200 μl	1.5	0.6

The error limits mentioned above are in accordance with the volume indicated on the instrument. These are obtained by using the instrument with distilled water at equilibrium, the ambient temperature of 20 °C while operating it smoothly and steadily.

Instructions for use:

1. Attach a compatible tip to the pipette tip cone.
2. Press the plunger all the way down. Dip 3-4 mm part of the tip in the given liquid and release the plunger. The liquid will fill up the tip.
3. Press the plunger all the way down again to dispense the liquid back in the reservoir.
4. While keeping the plunger pressed, touch the tip to the side of the reservoir to dispense any remaining liquid.
5. Release the plunger. Ignore the traces of liquid remaining in the tip.

The pipette is ready for use.

Repeat steps 2 to 4.

Designed and Manufactured by MICROLIT

629, Pakramau, Kursi Road,
Lucknow - 226026, UP, India

Call at: +91-8081461996, +91-9918625629

For queries: info@microlit.com

For feedback: customercare@microlit.com

www.microlit.com



EXPERIENCE PRECISION

MICROLIT LILPET PRO

In-built Tip Ejector

Two-step Plunger Movement

Fully Autoclavable

Fixed Volume



SPECIFICATIONS AND ORDERING INFORMATION	Model	Volume	Acc. $\pm\%$	CV $\pm\%$
	LILP-5	5 μl	2.5	1.5
	LILP-10	10 μl	1.2	0.8
	LILP-20	20 μl	1	0.5
	LILP-25	25 μl	1	0.4
	LILP-50	50 μl	1	0.4
	LILP-100	100 μl	0.8	0.3
	LILP-200	200 μl	0.8	0.3

The error limits mentioned above are in accordance with the volume indicated on the instrument. These are obtained by using the instrument with distilled water at equilibrium, the ambient temperature of 20 °C while operating it smoothly and steadily.

Instructions for use:

1. Attach a compatible tip to the pipette tip cone.
2. Press the plunger all the way down. Dip 3-4 mm part of the tip in the given liquid and release the plunger. The liquid will fill up the tip.
3. Press the plunger all the way down again to dispense the liquid back in the reservoir.
4. While keeping the plunger pressed, touch the tip to the side of the reservoir to dispense any remaining liquid.
5. Release the plunger. Ignore the traces of liquid remaining in the tip.

The pipette is ready for use.

Repeat steps 2 to 4.

Designed and Manufactured by MICROLIT

629, Pakramau, Kursi Road,
Lucknow - 226026, UP, India

Call at: +91-8081461996, +91-9918625629

For queries: info@microlit.com

For feedback: customercare@microlit.com

www.microlit.com



EXPERIENCE PRECISION