



OUR REACH



YOUR SATISFACTION IS OUR PRIORITY

Microlit products follow GMP norms and are manufactured under comprehensive QA (Quality Assurance) and QC (Quality Control). However, if you are dissatisfied with the operation of any of our products, get in touch with us or call your nearest Microlit dealer for free replacement.

CONTACT US

Microlit India

629 Pakramau, Kursi Road, Lucknow 226026, India
Phone: +91 9918625629, Email: info@microlit.com
www.microlit.com

Microlit USA

33 Wood Avenue South, Suite 600, Iselin, NJ 08830, USA
Phone: +1 732 321 0852, Email: info-usa@microlit.com
www.microlit.us

Microlit LATAM

Business Place Ibirapuera, Alameda dos Jurupis,1005,
Indianópolis, São Paulo, Brasil
Phone: +55 11 93269-5932, Email: info-latam@microlit.com

MICROLIT FAVEO | Carousel Stand for 6 Micropipettes




Microlit Faveo- Carousel Stand for six micropipettes is a cost effective blend of features and functionality.

Built by our in-house team of product design engineers, Microlit Faveo can hold both single channel and multichannel micropipettes.

The sleek design of the stand allows the pipettes to be conveniently stored in style. The rotating nature of this carousel stand allows the user to place or remove the pipettes easily from the stand. This helps the user in accessing the desired pipette with much ease.

The stable base of the stand assists in single handed operation. The sturdy design ensures safe keeping of the pipettes in a quite effective manner.


Made of polypropylene, this stand can be easily cleaned with soap and water, or common laboratory disinfectants.



Stand for 6 Micropipettes
Microlit Faveo is designed to hold up to 6 single channel micropipettes.




Can Store Both Single Channel and Multi Channel Micropipettes
This stand can hold up to 3 multichannel micropipettes, 6 single channel micropipettes or a combination of 3 single channel and 3 multichannel micropipettes.




Portable & Light Weighted
Due to its light weight and small size, this carousel stand can be easily placed in a laboratory without occupying much space.


360° Rotating Head
The rotating nature of Microlit Faveo allows the user to place or remove the pipettes easily from the stand.



Stable Base
The stable base of this carousel stand assists in single handed operation.



Easy to Clean
Microlit Faveo can be cleaned with soap and water, or common laboratory disinfectants.



Specifications & Ordering Information:

Model No. : MCS-6
Dimensions : Height: 323 mm, Base Dia: 195 mm
Material : Acrylonitrile butadiene styrene (ABS)
Autoclavability : Not Autoclavable
Color : White and Blue