

### ANACON LABORATORIES PVT. LTD.

ISO 9001:2015, ISO 14001:2015, OHSAS 18001:2007 Certified Organization  
 Recognized By Ministry of Environment & Forests (MoEF), New Delhi  
 Accredited By Quality Council of India by - NABET

Approved by : ♦ Food Safety & Standards Authority of India Under FSS Act ♦ Bureau of Indian Standards (BIS)  
 ♦ Agricultural & Process Food Products Export Development Authority (APEDA) ♦ Export inspection council (EIC) ♦ Agmark

Test Report No. : ALPL/17032020/03-1

dated 17.03.2020

Page 1 of 1

|  |   |                                    |
|--|---|------------------------------------|
| <b>Issued To :</b><br>M/s Microlit<br>629, Pakramau, Kursi Road,<br>Lukhnow (UP)-226026<br><b>Kind Attention :</b> Ms. Vidushi Srivastava<br><b>Contact No. :</b> 7073342594 | <b>Sample Inward No.</b> 2020/HN-335/1-1              | <b>Analysis Start</b> 06.03.2020   |
|  | <b>Inward Date</b> 06.03.2020                         | <b>Analysis End</b> 16.03.2020     |
|  | <b>Reference</b> ALPL/2019-20/Proposal/<br>G-3231(R1) | <b>Sample Category</b> General     |
|  | <b>Reference Date</b> 20.02.2020                      |                                    |
| <b>Sample Name</b><br>Bottle Top Dispenser   | <b>Sample Condition</b><br>Satisfactory               | <b>Quantity Received</b><br>1piece |
| <b>Sample Collected By</b><br>M/s Microlit   |   |                                    |


### TEST RESULTS

| S.N.      | Test Parameter  | Test Result             |
|-----------|---|-------------------------|
| <b>A.</b> | <b>Set 1</b>  |                         |
| 1.        | Visual inspection before testing  | No deformation observed |
| 2.        | Visual inspection after exposure to HF for 4 days at room temperature by filling up with HF | No change observed      |
| <b>B.</b> | <b>Set 2</b>  |                         |
| 1.        | Visual inspection before testing  | No deformation observed |
| 2.        | Visual inspection after exposure to HF for 4 days at room temperature by immersing in HF    | No change observed      |


**NOTES:** ● Please see watermark "Original Test Report" to confirm the authenticity of this report. ● Results shall be referred to tested sample(s) and applicable to tested parameters only.  
 ● Test report shall not be reproduced except in full without prior written approval of Anacon Labs. ● Liability of Anacon Labs is limited to invoiced amount only. ● Non-perishable and perishable sample(s) shall be disposed off after 30 days and 15 days respectively from the date of issue of Test Report, unless specified otherwise. ● Photograph attached.

**REMARKS :** As requested by the client, sample was tested for above parameters only. Above test are carried out in 50% concentrated hydrofluoric acid.

Verified By

  
 Swati Shrivastava  
 Technical Manager

Authorized Signatory

  
 ✓ Chinmay Garway  
 Director