



Micropipette Starter Kit

Add precision & convenience to
your pipetting needs



EXPERIENCE PRECISION



Product Code : MSK-04

Microlit Starter Kit Contains:

A. Set of Four RBO Series Micropipettes| Single Channel, Variable Volume

- 0.5-10 μ l
- 10-100 μ l
- 100-1000 μ l
- 1-10 ml

B. Tips & Tip Boxes

- Tip Box consisting 96 tips | Capacity upto 1000 μ l
- Tip Box consisting 96 tips | Capacity upto 100 μ l
- Tip Box consisting 96 tips | Capacity upto 10 μ l
- Packet of 10000 μ l capacity tips | Qty. 50 tips

C.Faveo | Carousel Stand for 6 Micropipettes

- Sleek carousel stand to provide convenient storage for micropipettes

Each micropipette is fully autoclavable and is individually calibrated in Microlit's NABL accredited lab certified as per ISO 17025 and is supplied with a calibration certificate conforming to ISO 8655 standards.

CONTACT US

★ Microlit INDIA

629 Pakramau, Kursi Road, Lucknow 226026, India
Phone: +91 9918625629, Email: info@microlit.com
www.microlit.com

★ Microlit LATAM

Business Place Ibirapuera, 9th Floor, Number 93,
Alameda dos Jurupis 1005 - Moema,
Indianópolis São Paulo, Brasil CEP - 04088-000,
Phone: +55 11 93269-5932, Email: info-latam@microlit.com

★ Microlit USA

33 Wood Avenue South, Suite 600, Iselin, NJ 08830, USA
Phone: +1 732 321 0852, Email: info-usa@microlit.com
www.microlit.us

Follow us

