


 MICROLIT

# MINI-STIR


DIGITAL MAGNETIC  
MINI STIRRER

- ONE TOUCH POWER ON/OFF BUTTON
- IP65 COMPLIANT SPILL PROOF DESIGN
- SPEED REGULATIONS UPTO 1500 RPM
- DIGITAL LED DISPLAY
- MOTOR LESS MAGNETIC STIRRER



# TECHNICAL SPECIFICATIONS

Specification	Motor less Magnetic Stirrer
Model Name	<b>MINI-STIR</b>
Model No.	<b>STIR-MI</b>
Stirring Position	<b>1</b>
Speed Increment / Decrement	<b>50 RPM</b>
Display	<b>Digital LED</b>
Working Plate Material	<b>304 Stainless Steel Grade</b>
Setting Dia	<b>105mm</b>
Maximum Stirring Capacity	<b>800ml H2O</b>
Recommended Stirrer Bar Size	<b>25x8mm PTFE Coated</b>
Pulse Mode	<b>30 seconds auto reverse direction</b>
Timer	<b>99 minutes and Infinity</b>
Last run Memory Save	<b>Yes</b>
Protection Class	<b>IP65 Compliant</b>
Permissible Ambient Temp	<b>5°C to 40°C</b>
Permissible Relative Humidity	<b>80%</b>
Voltage	<b>100-240 VAC</b>
Frequency	<b>50/60 Hz</b>
Weight	<b>0.3 Kg</b>
Dimension(WxDxH)	<b>162 x 12 x 110 mm</b>
Standard Accessories	<b>Stirrer Bar 1 Piece, 9V 1 Power Adapter, Operation Manual</b>



## CONTACT US

### Microlit INDIA

629 Pakramau, Kursi Road, Lucknow 226026, India  
Phone: +91 9918625629, Email: info@microlit.com  
[www.microlit.com](http://www.microlit.com)

### Microlit LATAM

Business Place Ibirapuera, 9th Floor, Number 93,  
Alameda dos Jurupis 1005 - Moema,  
Indianópolis São Paulo, Brasil CEP - 04088-000,  
Phone: +55 11 93269-5932, Email: info-latam@microlit.com

### Microlit USA

33 Wood Avenue South, Suite 600, Iselin, NJ 08830, USA  
Phone: +1 732 321 0852, Email: info-usa@microlit.com  
[www.microlit.us](http://www.microlit.us)

### Follow us

